Included in package:



1 Fabrication Dummy

Manufactured by:



7450 W. Lemhi Boise, ID 83709 USA (800) 819-5980 or (208) 429-0026 www.coyotedesign.com



CD103SFD

Fabrication Dummy

For CD103S, CD103SD Small Air-Lock and Small Deep Airlock



2-Year warranty against manufacturer defects, excessive wear or breakage

Made in U.S.A.

External Prosthetic Components



Advena Ltd

Pure Offices, Plato Close Tachbrook Park

Warwick, CV34 6WE, UK

COM CD103FD.revA.03102015





1 Fabrication Dummy



CD103SFD

Fabrication Dummy

For CD103S, CD103SD Small Air-Lock and Small Deep Airlock



2-Year warranty against manufacturer defects, excessive wear or breakage

Made in U.S.A.

External Prosthetic Components



Advena Ltd Pure Offices, Plato Glose Tachbrook Park Warwick, CV34 6WE, UK

CE CD103FD.revA.03102015

7450 W. Lemhi Boise, ID 83709 USA (800) 819-5980 or (208) 429-0026 www.coyotedesign.com

Manufactured by: **Couple**

Flexible Inner Socket with Coyote Designs Fabrication Dummy CD103SFD



For Small Air-Lock put fabrication CD103SFD dummy over anchor.



Pull preferred flexible plastic over socket and installed fabrication dummy using a vacuum system.



Pressing in around fabrication dummy will help remove any air pockets.



If needed use parachute cord or rope to pull plastic tight to dummies edge.



Let plastic cool before cutting fabrication dummy free.



Cut tooling free with cast saw, knife or utility knife.



Remove tooling from anchor.





Sand flexible inner down so anchor is almost 100% exposed

Pull a PVA bag over the anchor and cast. Glue your lock to anchor and start laminating.

Flexible Inner Socket with Coyote Designs Fabrication Dummy CD103SFD



For Small Air-Lock put fabrication CD103SFD dummy over anchor.



Pull preferred flexible plastic over socket and installed fabrication dummy using a vacuum system.



Pressing in around fabrication dummy will help remove any air pockets.



If needed use parachute cord or rope to pull plastic tight to dummies edge.



Let plastic cool before cutting fabrication dummy free.



Cut tooling free with cast saw, knife or utility knife.



Remove tooling from anchor.







Sand flexible inner down so anchor is almost 100% exposed



anchor and start laminating.