


PARTS IN THIS PACKAGE



- A. CQL Valve Base B. Valve Body C. Nail D. Felt Filter
E. Fabrication tool F. O-Ring

Manufactured By



419 N. Curtis Rd., Boise, ID 83706
(208) 429-0026
www.coyotedesign.com



CD105 CQL Valve

Weight limit 265 lbs



MADE IN USA

External Prosthetic Components



CD105.revA.02092016

Advena Ltd
Pure Offices, Plato Close
Tachbrook Park
Warwick, CV34 6WE, UK

RETRO-FITTING CQL™ VALVE



1. Mark the valve body location as close to the distal end of the socket as possible.

2. Typically valve is posterior



3. Use an 1/8" drill bit for pilot hole. Have flexible insert in socket when drilling pilot hole.



4. Remove insert before drilling 1/2" hole.



5. Drill a 1/2" for the valve body hole.



6. Carefully clean edges after drilling.



7. Ensure the Valve Base fits flush with the inside of the socket.



8. Run a small bead of Coyote Quick Adhesive around the edge of the valve base.



9. Carefully place the valve base in drilled hole. Hold base firmly in place.

10. Clean any excess adhesive.

11. Place O-ring on valve body and thread valve body into base.

12. Once adhesive is set place felt filter in bottom of the base.

13. Use a lock wrench to ensure the valve body is fully installed. Valve body should be hand tight.

If you have more questions please call
Coyote at (208) 429-0026 or toll free at (800) 819-5980.

LAMINATING THE CQL™ VALVE USING DUMMY



14. Mark valve location on inner PVA bag and cover mark with tape to prevent tearing.

15. Using provided nail attach lamination tool to mold.

16. Lay up and laminate as desired.



17. Carefully expose lamination tool and remove lamination tool.



18. Apply a bead of adhesive to the lip of the valve base and carefully place in opening of socket. Press

19. Apply a bead of adhesive around seam between valve base and socket.



20. Place black felt circle in valve base.



21. Carefully thread valve body into valve base. Using lock wrench tighten to hand tight.

The CQL Valve is an auto expulsion valve only

Results can only be guaranteed if Coyote Quick Adhesive is used.